

Transport Bottle

(#401000, 411000)

1. Description

SPL Transport Bottle for safe transport of pathogenic organisms and clinical specimens tightly sealed to prevent possible leakage.

2. Features

- ❖ Packaging for infectious substances (UN3373)
- ❖ Category A: An infectious substance which is transported in a form that, when exposure to it occurs, is capable of causing permanent disability, life-threatening or fatal disease in otherwise healthy humans or animals.
- ❖ Category B: An infectious substance which does not meet the criteria for inclusion in Category A.

3. Product Range

Cat. No.	Color (Cap/Body)	Material (Cap/Body)	External Dimension (mm)	Internal Dimension (mm)	Total Volume (ml)	Empty Weight (g)	Packaging
401000	Red / Natural	PP / PP	164 x 164 x 168	132 x 155	1,000	260	1 / 10
411000	Red / Natural	PP / PP	112 x 110 x 190	107.5 x 176.3	1,000	147	1 / 10

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs




7. Specifications & Technical Drawing

(1) Transport Bottle (#401000)

- Specifications

Dimensions				Material (Resin)	Packaging
Cap		Body			
Length	29.00 mm	Length	109.30 mm	▪Body: Polypropylene(PP) ▪Cap: Polypropylene(PP)	▪ Pieces/Sleeve: 1 ▪ Pieces/Case: 10
Width:	132.20 mm	Width:	151.00 mm		


- Technical Drawing



(2) Transport Bottle (# 411000)
 - Specifications

Dimensions		Material (Resin)	Packaging
Cap	Body		
Length 25.00 mm Width: 107.53 mm	Length 176.30 mm Width: 90.17 mm	▪Body: Polypropylene(PP) ▪Cap: Polypropylene(PP)	▪ Pieces/Sleeve: 1 ▪ Pieces/Case: 10

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
 To place an order, contact your local distributor or
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-TSPBO

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
 © SPL Life Sciences Co., Ltd., Version 2, May 23th, 2017

