

Insect Breeding Box

(#310075/310076/310077/310074)

1. Description

Insect Breeding boxes are ideal for breeding and observing micro-insects. Ventilation hole on the cover enables excellent air circulation and easy humidity control.

2. Features

- ❖ Easy to breed and observe micro-insects and their natural enemies
- ❖ Easy air circulation and humidity control
- ❖ Useful for room tests such as room breeding and drug screening
- ❖ Mesh type: stainless or nylon (please inquire)
- ❖ All products can be assembled with Incu Tissue (Cat. No. 310070)
- ❖ Optimal Frame (Cat. No. 310074)

3. Product Range

Cat. No.	Material	External Dimension w x l x h (mm)	Internal Dimension w x l x h (mm)	Ventilation Hole size (mm)	Mesh Pore Size (mm)	Hole Position	Packaging
310075	PS	72 x 72 x 100	65.4 x 65.4 x 98.2	40	0.053	Cap	4 / 120
310076	PS	72 x 72 x 100	65.4 x 65.4 x 98.2	40	0.053	Bottom	4 / 120
310077	PS	72 x 72 x 100	65.4 x 65.4 x 98.2	40	0.053	Slide&Bottom (3 Holes)	4 / 120
310074	PS	80.8 x 80.8 x 21.0	63.8 x 63.8 x 21.0	-	-	-	30 / 60

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs



SPL-TDS-IBB

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, July 21th, 2016

6. Specifications & Technical Drawing

- Specifications

Dimensions	Material (Resin)	Packaging
Hole: 40.00 Ø Length _{body bottom} : 65.00 mm Width _{body bottom} : 65.00 mm Length _{cap bottom} : 72.00 mm Width _{cap bottom} : 72.00 mm Height: 99.20 mm	Crystal-grade polystyrene (PS)	▪ Pieces/Sleeve: 4 ▪ Pieces/Case: 120

- Technical Drawing

Cat. No. 310075 (Cap)	Cat. No. 310076 (Bottom)	Cat. No. 310077 (3T)

www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-IBB

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, July 21st, 2016