

Bacteria Culture Tube

(#59015/59050)

1. Description

During cultivation of pathogenic microorganisms, researchers can handle the tube while protecting the sample from infectious contaminations.

2. Features

- ❖ Sampling of the culture supernatant can be performed without infection during the incubation of pathogenic microbial cells such as mycobacterium tuberculosis, etc.
- ❖ Syringe filter is mounted on the cap to minimize contamination from pathogenic microorganisms.
- ❖ 15/50 ml conical tube size.
- ❖ Suitable for shaking incubation.

3. Product Range

Cat.No.	Material (Tube/Cap/Filter)	Cap Type	External Dimension d x h (mm)	Total Volume (ml)	Sterile	Packaging
59015	PP/HDPE/PVDF	Filter	17.00 x 140.50	15.00	+	5 / 50
59050	PP/HDPE/MCE	Filter	30.00 x 137.90	50.00	+	5 / 50

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Specifications & Technical Drawing

(1) Bacteria Culture Tube (#59015)

- Specifications

Dimensions		Material (Resin)	Packaging
Cap	Diameter: 21.30 mm Height: 33.00 mm	High-density polyethylene (HDPE)	▪ Pieces/Sleeve: 5 ▪ Pieces/Case: 50
Body	Diameter: 16.90 mm Height: 119.20 mm Volume: 15.00 ml	Polypropylene(PP)	

- Technical Drawing




(2) Bacteria Culture Tube (#59050)

- Specifications

Dimensions		Material (Resin)	Packaging
Cap	Diameter: 34.04 mm Height: 36.40 mm	High-density polyethylene (HDPE)	▪ Pieces/Sleeve: 5 ▪ Pieces/Case: 50
Body	Diameter: 29.30 mm Height: 113.90 mm Volume: 50.00 ml	Polypropylene(PP)	

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr