

Tube Top Filter Unit (#51850/51851)

1. Description

The SPL Tube Top Filter Unit can directly filter reagents contained in a 50 ml Conical Tube, and since the filtered reagents are immediately placed in another 50 ml Conical Tube, which provides the convenience of eliminating the need for a secondary step of transferring the sample.

It is very useful for sterile-filtration of media or reagent at a small volume, and for separation of cells.

2. Features

- ❖ Individual Packing
- ❖ 2 pore size: 0.22, 0.45 μm

3. Functional Test

Non-cytotoxic : Passed

4. Product Range

Cat. No.	Material (Tube/Cap/Membrane)	Pore size (μm)	Tube Volume (ml)	Sterile	Packaging
51850	PP/PS/PES	0.22	50.00	+	1 / 25
51851	PP/PS/PES	0.45	50.00	+	1 / 25

5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2008 Quality Management System (K-QA-Q031478) in Korea.

6. Product Photographs




7. Specifications & Technical Drawing

- Specifications

Dimensions		Material (Resin)	Packaging
Height	138.20 mm	Polystyrene (PS)	▪ Pieces/Sleeve: 1 ▪ Pieces/Case: 25
Bottom Diameter	27.89 mm	Polypropylene (PP)	
		High Density Polyethylene(HDPE) Polyethylsulfone(PES)	

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
To place an order, contact your local distributor or
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-GridSticker

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 1, June 15th, 2017

